



SURFACE SANITISING WIPES

Suitable for a wide range of critical cleanroom, laboratory and manufacturing environments.



LabService Providers

Douglas Guilfoyle, Sales Director

087 0669949

doug@labserviceproviders.ie



LabService
Providers

Surface Disinfectant Wipe

- Quat and PHMB free.
- Excellent efficacy & kill rate for both gram positive & gram negative organisms.
- Highly effective broad spectrum efficacy proven by independent testing.
- Class leading efficacy at lower temperatures - 5 - 10 degrees celsius - and under 'dirty' conditions.
- Registered under Annex 1 of the Biocides Products Regulations (98/8/EEC) under 3 Product Types: PT 2 Private and Public Healthcare environments, PT 3 Veterinary Hygiene and PT 4 Food and Feed area disinfectants.
- Dept Agriculture approved PCS no.
- Low corrosion on metal and steel surfaces due to low salt content.
- Non oxidising - optimum continuous efficacy in environments where organic matter is present on surfaces.
- Sheet size: 200 x 180 mm.
- 1 Litre bucket / 1500 sheets.



70/30 IPA Wipe

- 70% IPA disposable surface wet wipe designed for the effective disinfection of all surfaces,
- Prevents cross contamination in all healthcare areas.
- Tested and approved to European norms EN1276/EN14348/EN13727 for bactericidal efficacy in medical environments.
- Effective kill rates against all potentially lethal micro-organisms including MRSA, HIV, C-Difficile, Enterococcus Hirae.
- Tested and conform to EN13624 fungicidal activity,
- Also ideal for organic food production as the IPA wipe leaves no residue.
- White aquaspun polypropylene substrate (23gsm),
- Sheet Size: 200mm x 200mm.
- 200 Wipes Per Tub.

